


Summary of Results from  
Group 4  
Lactase Downstream Processing

Evreux 23. – 27.03.2026

# Group 4 + presentation of group members

Manon Sala  
French



Fabian Klingel  
German



Rasa Rudzevičienė  
Lithuanian



Hana Melcherová  
Slovakian




# Evaluation Steps of our Experiment

1. **PARTICLE FREE LYSATE – o-NPG-test (enzymatic activity assay)**
2. ANION EXCHANGE CHROMATOGRAPHY (AEXC) – o-NPG + photometer (enzymatic activity assay)
3. SDS – PAGE – gel + calibration line (purity, molecular weight – lactase)
4. RAW MILK TEST – colour change (conc. of glucose)
5. EVALUATION

# O-NPG-TEST of 1:10 diluted Particle Free Lysate

Test 50  $\mu\text{L}$  supernatant + 500  $\mu\text{L}$  o-NPG-solution.

- **Observation:** Yellow
- **Evaluation:**  $\beta(1\rightarrow4)$  bond was broken down
- **Conclusion:** lactase is present.  
→ We can start AEXC.



# Evaluation Steps of our Experiment


1. PARTICLE FREE LYSATE – o-NPG-test (enzymatic activity assay)
2. **ANION EXCHANGE CHROMATOGRAPHY (AEXC) – o-NPG + photometer (enzymatic activity assay)**
3. SDS – PAGE – gel + calibration line (purity, molecular weight – lactase)
4. RAW MILK TEST – colour change (conc. of glucose)
5. EVALUATION

# ANION EXCHANGE CHROMATOGRAPHY

O-NPG-test measured by photometer – enzymatic activity measurement

Test 50  $\mu\text{L}$  eluate + 800  $\mu\text{L}$  o-NPG-solution.

- **Observation:** one sample turned Yellow
- **Evaluation:** only one of the samples contains lactase
- **Conclusion:** the highest lactase activity was in eluate 25%



# Evaluation Steps of our Experiment

1. PARTICLE FREE LYSATE – o-NPG-test (enzymatic activity assay)
2. ANION EXCHANGE CHROMATOGRAPHY (AEXC) – o-NPG + photometer (enzymatic activity assay)
3. **SDS – PAGE – gel + calibration line (purity, molecular weight – lactase)**
4. RAW MILK TEST – colour change (conc. of glucose)
5. EVALUATION

# SDS - PAGE

## MARKER

**S1:** lysate (dilution factor 10)

**S2:** loading flow through


**S3:** washing step

**S4:** eluate 10%

**S5:** eluate 25%


**S6:** eluate 100%

- Conclusion: by AEXC we reduced the number of proteins in our sample to an uncountable amount of proteins (starting with 3000)




# SDS - PAGE

- We identified the lactase subunit (116kDa) in our gel using marker
- Lactase Band is a little bit above 100 marker



# Evaluation of the SDS – PAGE

- Identification of molecular weight of lactase subunit with your calibration line
- **Migration rate: 5,6cm**
- **Molecular weight: 117kDa**




# Evaluation Steps of our Experiment

1. PARTICLE FREE LYSATE – o-NPG-test (enzymatic activity assay)
2. ANION EXCHANGE CHROMATOGRAPHY (AEXC) – o-NPG + photometer (enzymatic activity assay)
3. SDS – PAGE – gel + calibration line (purity, molecular weight – lactase)
4. **RAW MILK TEST – colour change (conc. of glucose)**
5. EVALUATION

# RAW MILK TEST

- **Observation:** barely any colour change after 20 and 40 mins. After 70 mins and applying heat, sample turned green
- **Evaluation:** glucose level of approximately 2.8 mmol/L
- **Conclusion:** to create lactose free milk we would have to change the conditions so the enzyme can work better.



# Evaluation Steps of our Experiment

1. PARTICLE FREE LYSATE – o-NPG-test (enzymatic activity assay)
2. ANION EXCHANGE CHROMATOGRAPHY (AEXC) – o-NPG + photometer (enzymatic activity assay)
3. SDS – PAGE – gel + calibration line (purity, molecular weight – lactase)
4. RAW MILK TEST – colour change (conc. of glucose)
5. **EVALUATION**

# EVALUATION

- 1. TASK:** SDS-PAGE and AEXC both showed the same results. 25% eluate had the highest enzymatic activity and showed the lactase subunit band.
- 2. TASK:** We would use the 25% eluate, we would have to further purify it by chromatography.

Thank you for your  
attention!